A little different than normal...

We will do the first part of the lecture in our desks, and then the rest of this PowerPoint will be done as you do the lab activity.

N50 – Titrations Our last lecture of the year!

Target:

I can set up and perform a titration.

What is titration?

A way to determine the concentration of an unknown substance.

- Uses the fact that acids and bases react with each other in "neutralization reactions"
- At the point where the neutralization reaction is finished # moles Acid = # moles Base

Key Terms

Titrand

The unknown solution you are interested in

Titrant

The solution with the known concentration

Equivalence Point

The point at which all the titrand has reacted with the titrant. # Moles Acid = # Moles Base

End Point


The point at which your titration <u>seems</u> finished during the lab – a color change happens for example

How do you know you reached the end point?


Use an **INDICATOR**

Turns colors based on pH – can show you visually when you have reached the end point.


Pick the right indicator!




Pick the right indicator!




Lab Set Up





Careful! Don't go too fast!



So... Known [base] & unknown [acid]

Volume
Base used
to get to
end point
From

burette

during lab

Convert
Vol → Moles
using M = mol/L

Moles
Base used
to get to
end point

At end point
(assuming you did a
good job titrating!)

Moles Acid = Moles Base

Acid in unknown solution

Moles

Convert
Moles → Molarity
using M = mol/L

Concentration of unknown Acid

So... Known [acid] & unknown [base]

Volume
Acid used
to get to
end point
From

burette

during lab

Convert
Vol → Moles
using M = mol/L

Moles
Acid used
to get to
end point

At end point (assuming you did a good job titrating!)

Moles Acid = Moles Base

Convert
Moles → Molarity
using M = mol/L

Concentration of unknown Base

Moles

Base in

unknown

solution

Some things to be careful of...

mL versus L

Stoichiometry - is it a 1:1 ratio H⁺: OH⁻?
Or is it 1:2, or 2:1, or 2:3, etc
1 mol NaOH = 1 mol OH⁻
1 mole Ca(OH)₂ = 2 mol OH⁻

 End point and equivalence point are only identical if your titration is absolutely perfect. It never is, there are lab errors!

Lab Activity Portion of Lecture

Instead of a practice problem just on paper, we are going to have our practice problem be an actual titration! On WS #14

Question #1: What is our titrand?

Question #2: What is our titrant?

```
Titrand = HCI – unknown []
Titrant = NaOH – known [], 0.10 M
```


Lab Set Up

- 1) Make sure everything is rinsed with distilled water, and then rinse the burette with the titrant (NaOH).
- 2) Clamp burette into burette clamp onto a ring stand
- 3) Fill burette with NaOH with a known concentration (0.10 M)
- 4) Put a beaker under burette and slowly open valve, letting some NaOH out until bottom of meniscus is reading at an easy to read value. Careful! Make sure that the entire tip of the burette is filled with NaOH.


Lab Set Up continued...

- 5) Using a graduated cylinder, put a known volume of your titrand (unknown concentration HCL) into an Erlenmeyer flask.
- 6) Add a small amount of appropriate indicator to the flask (phenolphthalein).

Question #3: Which indicator should we pick? Our unknown will be in the 7-10 range.



Question #4: What color shift do you expect to see for this indicator? From _____.



How many trials to do?

Do FOUR trials (typically)

1st Trial – Rough trial – "Quick and dirty"

Just a rough estimate so you have an idea of when you need to start slowing down. DON'T include this trial when averaging your data!

2nd – 4th Trials – Real ones – Be careful!

Set up your Data Table

| | Titration of Unknown HCl Solution with Phenolphthalein | | | | | |
|------------------------------------|--|----------|--|---------------------|----------|--|
| | [] of Titrant: Volu | | | me of Titrand Used: | | |
| | Rough Trial | Trial #1 | | Trial #2 | Trial #3 | |
| Name | | | | | | |
| Burette Starting Volume (mL) | | | | | | |
| Burette Ending Volume (mL) | | | | | | |
| Volume of Titrant Used (mL) | | | | | | |

Every group needs to do one rough trial, and each person will do a real trial. Four people in a group = Four trials. Make the right number of columns in your table!

Set up your Data Table

| | Titration of Unknown HCl Solution with Phenolphthalein | | | | | |
|--------------|--|----------|-------------------------------|----------|----------|--|
| | [] of Titrant: 0.10 M | | Volume of Titrand Used: 5 mL | | | |
| | Rough Trial | Trial #1 | | Trial #2 | Trial #3 | |
| Name | | | | | | |
| Burette | | | | | | |
| Starting | | | | | | |
| Volume (mL) | | | | | | |
| Burette | | | | | | |
| Ending | | | | | | |
| Volume (mL) | | | | | | |
| Volume of | | | | | _ | |
| Titrant Used | | | | | | |
| (mL) | | | | | | |

Titration Lecture Videos

- 1) What is a titration? FuseSchool https://youtu.be/tlbD8MG1qMM
- 2) Setting up and Performing a Titration. CarolinaBiological https://youtu.be/sFpFCPTDv2w
- 3) How to Prepare a Burette for a Titration. Wits University https://youtu.be/Lr1nLTCqZvM
- 4) How to Read the Volume off a Burette. Wits University https://youtu.be/qdmp4_Nwd-Q
- 5) What is a Titration and how is it performed? Wits University https://www.youtube.com/watch?v=YqfvRBJ-iPg
- 6) Acid Base Equilibrium. Bozeman Science https://youtu.be/l5fk7HPmo5g

Set up your lab station

Burette's are already clamped in for you, and filled with NaOH


- not super safe to fill them, so I did it for you.

Perform Rough Trial

| | Titration of Unknown HCl Solution with Phenolphthalein | | | | | |
|------------------------------------|--|-------|-------------------------------|----------|----------|--|
| | [] of Titrant: 0.10 M | | Volume of Titrand Used: 5 mL | | | |
| | Rough Trial | Trial | #1 | Trial #2 | Trial #3 | |
| Burette Starting Volume (mL) | | | | | | |
| Burette Ending Volume (mL) | | | | | | |
| Volume of Titrant Used (mL) | | | | | | |

Each person does their trial

| | Titration of Unknown HCl Solution with Phenolphthalein | | | | | |
|------------------------------------|--|----------|-------------------------------|----------|----------|--|
| | [] of Titrant: 0.10 M | | Volume of Titrand Used: 5 mL | | | |
| | Rough Trial | Trial #1 | | Trial #2 | Trial #3 | |
| Burette Starting Volume (mL) | | | | | | |
| Burette Ending Volume (mL) | | | | | | |
| Volume of Titrant Used (mL) | | | | | | |








When your titration turns bright pink, and the professor starts walking towards you...







The perfect game doesn't exi-

Calculate how many moles of NaOH you used

Molarity = Moles / Liters From burette in the lab

Moles NaOH used = Volume used x Molarity NaOH

Moles = L x mol L given to you

Moles NaOH = x mL 1L Y mol 1000 mL 1 L

Calculate the unknown concentration of the acid

At End Point -> Moles NaOH = Moles HCI

Molarity = Moles / Liter

Same as mol NaOH used!

Molarity Acid = Moles Acid Liters Acid Used

The amount in your Erlenmeyer flask!

Come check what your unknown concentration of Acid was!

Calculate the % error for each person's acid – let's see which person had the most accurate titration per group!

Average your group member's answers together – report % error of averaged data on the whiteboard – let's see which group had the best titration skills!

YouTube Link to Presentation

► https://youtu.be/6owm822vyhl