## Dougherty Valley HS Chemistry Types of a Reactions Lab



Name: Period: Seat#:

Required Sections: (Refer to R-5 for guidelines and requirements. Make note of any specific changes given by your teacher in class)

Prelab: All written in your lab notebook - Answer Pre-Lab Questions, Materials, Reagent Table, Procedures,

**During Lab:** Data Tables are part of the actual Lab Handout this time!

Post-lab: - The written equations on the actual Lab Handout are the post lab questions this time!

## Introduction

Most reactions fall into one of five categories that we have discussed in class. In this lab, you will be completing several different types of reactions. For EACH reaction you investigate, you must predict the products, identify the type of reaction, and write both the word equation AND the balanced formula equation.

## **Pre-Lab Questions**

- 1. What are the five main categories of reactions?
- 2. Describe how to predict the products for each type of reaction.

Reaction #1				
Materials  ■ CuSO₄  ■ Large paperclip ■ Test tube	<ol> <li>Procedure</li> <li>Add a little CuSO<sub>4</sub> to a test tube. CuSO<sub>4</sub> is toxic. Handle it with care!</li> <li>Open a paper clip (source of Fe, iron) and hang it over the edge of a test tube reaching ~ 2cm into the CuSO<sub>4</sub>. Observe the paper clip for 5 to 8 minutes.</li> <li>Remove the paper clip and place it on a piece of white scratch paper. When finished, throw out the paper clip and recycle the copper sulfate in the container at the front desk.</li> </ol>			
Observations				
Type of Reaction				
Word Equation				
Balanced Equation with phases included				
Net Ionic Equation <u>(if</u> <u>applicable)</u>				

Reaction #2				
Materials  ■ 0.15M SrCl₂  ■ 0.25M Na2CO₃  ■ One 100mL beaker  ■ Erlenmeyer flask  ■ Filter paper  ■ Funnel  ■ Wash bottle with distilled	Procedure  1) Using a graduated cylinder, measure 15ml of SrCl <sub>2</sub> and put into the 100mL beaker.  2) Measure 15mL of Na <sub>2</sub> CO <sub>3</sub> and add to the SrCl <sub>2</sub> in your 100mL beaker.  3) Record your observations in the chart below.			
Observations				
Type of Reaction				
Word Equation				
Balanced Equation with phases included				
Net Ionic Equation <u>(if</u> <u>applicable)</u>				

Reaction #3			
Materials     Steel wool     Bunsen burne     Tongs	Procedure  1) Burn a small sample of steel wool (source of iron) over the Bunsen burner.  r		
Observations			
Type of Reaction			
Word Equation (assume Fe <sup>3+</sup> is formed)			
Balanced Equation with phases included (assume Fe³+ is formed)			
Net Ionic Equation <u>(if</u> <u>applicable)</u>			

Reaction #4				
<u>Materials</u>		Procedure  1) Francisco the "Decemposition Reactions of Acide" on the gold never on your		
<ul> <li>Gold paper reading on lab table</li> </ul>		1) Examine the "Decomposition Reactions of Acids" on the gold paper on your table. Carbonic acid (H <sub>2</sub> CO <sub>3</sub> ) is a component of most carbonated sodas.		
Observations				
Type of Reaction				
Word Equation				
Balanced Equation with phases included				
Net Ionic Equation <u>(if</u> <u>applicable)</u>				